


# ESHERO

Beijing Eternity Electronic Technology Co., LTD.



## At3

### Anesthesia Machine

All the contents of this product such as model, configuration and discription are for reference only, please refer to real products base on your need.



# At3 Anesthesia Machine

- Integrated design to create a simple and comfortable surgical space;
- Electric drive, Electric control, no Driving gas required;
- Meet different scenarios, VCV, PCV, VS, manual ventilation mode;
- Intelligent ventilation mode fill in weight automatic update parameters, simple operation;
- 8 inch touch screen;
- The minimum VTl can be 5ml;
- Three waveforms can be displayed on the same screen;

